

MSHOT MSX2 VS Olympus DP27



MSX2	DP27
12.0MP	5.0MP
1 inch	2/3"
15fps	15fps
380~1000nm	~
100us~1s	50us~8s
8bit/12bit	12bit
USB3.0	USB3.0
MSHOT software	MSHOT software
Windows XP, 7,8 ,10	Windows XP,7,8
Much lower	Around \$ 12,000




MSX2 Scientific Camera Benefits

- 1 inch large size sensor offer largest filed of view.
 - Super high 12 megapixels.
 - Rapid 15fps frame rate for full resolution.
- Spectral response 380nm~1100nm offer high light sensitivity.
 - Low noise for pure image background.
 - Competitive price.

Olympus CX31 - 10X -Objective- Pathology Sample (CD5-HE-CTC Stain)

MSX2

DP27



MSX2

DP27

