

# MshOt

## ML11-II Features:

- Exquisite design, field of view is wide and clear.
- 3W LED, Pre-set lamp center, brightness continues adjustable.
- Entry level microscope



## ML11-II

## Specification

| Item                     | Specification   |
|--------------------------|---|
| Eyepiece                 | WF10X/18  |
| Head                     | 30°inclination trinocular, P.D. 50mm-75mm, single diopter adjustable, fix light division R:T:80:20<br>30°inclination binocular, P.D. 50mm - 75mm, single type diopter adjustable (optional) |
| Objectives               | Achromatic objective 4X/0.10, W.D.: 37.5mm<br>Achromatic objective 10X/0.25,W.D: 6.55mm<br>Achromatic objective 40X/0.65,W.D: 0.669mm<br>Achromatic objective 100X/1.25, W.D: 0.198mm       |
| Nosepiece                | Quadruple nosepiece with inward tilt  |
| Focusing                 | 155mmx142mm double layer mechanical, moving range 76mmX50mm , accuracy 0.1mm  |
| Stage                    | Low position hand wheel coaxial coarse and fine adjustment with limit stopper and lock device, stage height movement (coarse movement stroke: 25mm),fine focusing accuracy 0.002mm          |
| Transmitted illumination | 3W LED, Pre-set lamp center, brightness continues adjustable<br>Abbey condenser,N.A.1.25, with changeable aperture diaphragm<br>Filters: blue, yellow and green                             |
| Camera port              | 1XC   |
| Power line               | External set wide voltage power adapter, input 100V-240V-AC50/60HZ, outputDC12V2A   |

Guangzhou Micro-shot Optical Technology Co. Ltd

Add: Room 506 Vanke Cloud Building A / 1933 Huaguan Road Guangzhou / China

WEB:www.m-shot.com / sales@m-shot.com / TEL:0086-20-38262481



# MshOt

## ML11-II

## Biological Microscope

